



by 🔱 3D SYSTEMS

Biocompatible 3D Printing Materials

Company

Nextdent B.V. was founded in 2012 in The Netherlands as a subsidiary company under Vertex Global Holding. The aim was to complement the dental materials portfolio of its affiliate Vertex-Dental B.V., which has a 79-year track record in research, development, regulatory and lean production of biocompatible dental materials. In January 2017, Vertex Global Holding merged with 3D Systems, a company that provides comprehensive 3D products and services, including 3D printers specifically designed for dental applications. Since the merger the unsurpassed NextDent material product range is being developed, registered, produced and sold through Vertex-Dental B.V., Soesterberg, The Netherlands.

Our Mission

The combination of 3D Systems transforming Figure 4™ 3D printing technology and NextDent's revolutionary materials will enable us to lead innovation in digital dentistry. We will leverage the unique opportunity to integrate our advanced materials and automated print workflows into transformative digital production solutions delivered to the dental practice and laboratory. Users of these solutions are able to produce trays, models, drilling templates, dentures, orthodontic splints, crowns and bridges with enhanced speed, precision and efficiency and lower cost compared to conventional procedures.

Certified and Biocompatible 3D Printing Materials

It is the sheer choice you have in printable materials for a wide range of dental applications that makes NextDent materials unique. Our products are classified in accordance with the medical device directive 93/42/EEC, making them suitable for local registrations as well. We offer material in a spectrum of colors, depending on the daily use and the preferences of our customers.

Experience You Can Count On

3D Systems has leveraged its 30 years of 3D printing experience in for the development of a number of dental 3D printers that offer high precision and detail resolution. They are used for the processing of different dental materials, from metal (CoCrMo alloy) to biocompatible resins. Our dental solutions are designed for use in dental laboratories, making production methods faster, easier and more effective.





Denture 3D+

NextDent Denture 3D+ is a biocompatible Class IIa material suitable for printing all types of removable denture bases. This material has significant lower shrinkage compared to standard PMMA denture base materials. Shrinkage can be compensated by using dental software resulting in excellent fitting denture bases. This material has excellent mechanical properties and is comparable to conventional denture base materials. Available in 5 new colors: Dark Pink, Light Pink, Opaque Pink, Red Pink and Translucent Pink.

Property	Requirement	Result	ISO standard
Ultimate flexural strength	≥65MPa	84	ISO 20795-1
Flexural modulus	≥2000 MPa	2383	ISO 20795-1
Sorption	≤32 µg/mm³	28	ISO 20795-1
Solubility	≤1,6 µg/mm³	0,1	ISO 20795-1
Residual monomer	≤2,2%(w/w)	<0,1	ISO 20795-1
Biocompatibility	Non-cytotoxic	Comply	
	Non-mutagenic	Comply	
	Not induce any erythema or edema reactions	Comply	ISO 10993-1
	Notasensitizer	Comply	
	Not cause systemic toxicity	Comply	





NextDent C&B Micro Filled Hybrid is a biocompatible Class IIa material developed for crowns and bridges and usable for denture teeth. The balance between inorganic fillers and the resin give the material its high strength and wear resistance. The material is easy to finish and polish, and can be stained with all types of composite staining kits. Due to the perfect balance between opacity and translucency the printed crown blends in perfectly between the existing teeth. Available in colors BL, N1, N1.5, N2, N2.5, N3 and T1.

Property	Requirement	Result	ISO standard
Flexural strength	≥50MPa	107	ISO 10477
Sorption	≤60µg/mm³	54	ISO 10477
Solubility	≤12,5 µg/mm³	5,9	ISO 10477
Biocompatibility	Non-cytotoxic	Comply	ISO 10993-1
	Non-mutagenic	Comply	
	Not induce any erythema or edema reactions	Comply	
	Notasensitizer	Comply	
	Not cause systemic toxicity	Comply	



Try-In

NextDent Try-In is a biocompatible Class I material suitable for printing try-in devices, a baseplate combined with the individual designed tooth setup, to check bite registration and occlusion. Available in colors TI0, TI1, TI2.


Property	Requirement	Result	ISO standard
Flexural modulus	≥1500 MPa	1882	ISO 20795-1
Biocompatibility	Non-cytotoxic	Comply	
	Non-mutagenic	Comply	ISO 10993-1
	Notasensitizer	Comply	



Tray

NextDent Tray is a biocompatible Class I material designed to print individual impression trays. The material distinguishes itself through high printing speed and accuracy. NextDent Tray resin enables technicians to make even the most complex trays in a matter of minutes. The printed trays are rigid and suitable for all types of impression material, making high-quality, high-precision impressions possible. Available in colors Blue and Pink.

Property	Requirement	Result	ISO standard
Ultimate Flexural strength	≥50 MPa	81	ISO 20795-1
Flexural modulus	≥1500 MPa	2015	ISO 20795-1
Residual monomer	≤2,2%(w/w)	<0,1	ISO 20795-1
Biocompatibility	Non-cytotoxic	Comply	
	Non-mutagenic	Comply	ISO 10993-1
	Notasensitizer	Comply	



SG

NextDent SG is a biocompatible Class I material, developed for the printing of Surgical Guides for implant surgery use. Because of the high precision of this material it is easy to insert drill sleeves, directly after printing. Enabling even greater precision during surgery. The properties of NextDent SG will not change by using disinfectants. In addition, the material can also be sterilized using standard autoclave protocols. The use of an autoclave does not affect the dimensional stability, therefore NextDent SG can be used in every operating theatre. Available in color Translucent Orange.

Property	Requirement	Result	ISO standard
Ultimate Flexural strength	≥50 MPa	85	ISO 20795-1
Flexural modulus	≥1500 MPa	2118	ISO 20795-1
Residual monomer	≤2,2%(w/w)	<0,1	ISO 20795-1
Biocompatibility	Non-cytotoxic	Comply	
	Non-mutagenic	Comply	ISO 10993-1
	Notasensitizer	Comply	



Model 2.0

NextDent Model 2.0 is characterized by its high degree of accuracy, making this material suitable for detailed master prosthodontic and orthodontic models where high precision is needed. The models show highly visual details due to color and opacity and has an ideal surface for scan impressions. The accurate printed models are the perfect base for creating your dental work piece. Available in colors Peach, White and Grey.

Property	Requirement	Result	ISO standard
Hardness shore D	≥80 Shore D	84	ISO 178



Gingiva Mask

NextDent Gingiva Mask is a flexible material that can be used in combination with the model material. This makes it possible to print parts of the model that need a certain flexibility, such as Gingiva Masks on implant models. Available in color Pink.

Property	Requirement	Result	ISO standard
Shore A hardness	60-75	68	ISO 10139-2
Elongation@ break	40-60%	153	ISO 527-1 ISO 527-2



Model Ortho

NextDent Model Ortho is a material suitable for printing models specifically used in 'Vacuum Molding' applications. This material is faster and easier to print compared to other model materials. This material can be printed faster compared to Model 2.0. Available in color Opaque Beige.

Property	Requirement	Result	ISO standard
Hardness shore D	≥80 Shore D	82	ISO 178



Ortho Clear

NextDent Ortho Clear is a is a biocompatible Class IIa material for all types of splints and retainers. It is a clear and thereby aesthetically-pleasing 3D print material, characterised by high strength and fracture thoughness. Available in color Clear.

Property	Requirement	Result	ISO standard
Ultimate flexural strength	≥50 MPa	70	ISO 20795-2
Flexural modulus	≥1300 MPa	1596	ISO 20795-2
Max. stress intensity factor	≥ 1,1 MPa m ^{1/2}	1,6	ISO 20795-2
Total fracture work	≥250J/m²	675	ISO 20795-2
Sorption	≤65µg/mm³	58	ISO 20795-2
Solubility	≤5,0µg/mm³	4,6	ISO 20795-2
Residual monomer	≤5,0% (w/w)	<0,1	ISO 20795-2
Biocompatibility	Non-cytotoxic	Comply	ISO 10993-1
	Non-mutagenic	Comply	
	Not induce any erythema or edema reactions	Comply	
	Notasensitizer	Comply	
	Not cause systemic toxicity	Comply	



Ortho Rigid

NextDent Ortho Rigid is a biocompatible Class IIa material developed for digital manufacturing of splints. In combination with suitable software, it is possible to easily design and print splints. Available in Transparent Blue.

Property	Requirement	Result	ISO standard
Ultimate flexural strength	≥50 MPa	78	ISO 20795-2
Flexural modulus	≥1500 MPa	2075	ISO 20795-2
Max. stress intensity factor	≥1,1 MPa m ^{1/2}	1,1	ISO 20795-2
Total fracture work	≥250J/m²	262	ISO 20795-2
Sorption	≤32 µg/mm³	20	ISO 20795-2
Solubility	≤5,0 µg/mm³	0,8	ISO 20795-2
Residual monomer	≤5,0 % (w/w)	<0,1	ISO 20795-2
	Non-cytotoxic	Comply	
	Non-mutagenic	Comply	
Biocompatibility	Not induce any erythema or edema reactions	Comply	ISO 10993-1
	Notasensitizer	Comply	
	Not cause systemic toxicity	Comply	



Ortho IBT

NextDent Indirect Bonding Tray is a biocompatible Class I material for orthodontic applications. Making use of the right dental software, you can plan the exact location of orthodontic brackets and design the Indirect Bonding Tray. Due to the flexible characteristics of the printed indirect bonding trays, orthodontist can easily place all the brackets at once, saving time at the chair. Available in color Clear.

Property	Requirement	Result	ISO standard
Shore A hardness	75-90	85	ISO 10139-2
Elongation@ break	10-20%	17	ISO 527-1 ISO 527-2
Biocompatibility	Non-cytotoxic	Comply	
	Non-mutagenic	Comply	ISO 10993-1
	Not a sensitizer	Comply	



Cast

NextDent Cast is a easy burn out 3D printing material, residue-free and suited for all kinds of purposes. Whether you need a support structure, a frame or orthodontic devices. Simply design and print making use of our castable material and the investment material that is recommended. Casting the parts in the metal you require. Available in color Purple.

Property	Requirement	Result	ISO standard
Flexural strength	≥60 MPa	85	ISO 178
Flexural modulus	≥1500 MPa	2193	ISO 178
Hardness shore D	≥80 Shore D	83	ISO 178



NextDent[™] 5100

The NexDent 5100 3D printer powered by revolutionary Figure 4[™] technology combined with NextDent's broad portfolio of dental materials addresses multiple indications, resulting in unparalleled speed, accuracy, repeatability, productivity, and total cost of operation.

The NextDent 5100 facilitates high-speed 3D printing for production of dental appliances and sacrificial castings. This revolutionary solution features an industry-defining value proposition that combines best-inclass speed and performance at a price point that is accessible to virtually all labs and clinics.

The range of dental materials, advanced print technology and compatibility with leading dental software is transforming dental workflows, enabling dental labs and clinics to produce trays, models, surgical guides, dentures, orthodontic splints, crowns and bridges with enhanced speed, precision and efficiency and lower cost.

Specification	Value
Related voltage	110V100-240VAC, 50-60HZ, 4.0A
Fuse	T5A/250V
Dimension 3D printer	W42.6xD48.9xH97.1cm
Dimension 3D printer with pedestal	W68.1 x D70.4 x H135.6 cm
Weight	35 Kg/54 Kg with pedestal



LC-3DPrint Box

The LC-3DPrint Box is a revolutionary UV light box, suitable for post-curing 3D printing materials. The new LC-3DPrint Box is equipped with 12 UV light bulbs strategically placed inside the box. This ensures that a product is illuminated from all sides, which results in a quick and uniform curing cycle. The spacious interior (Ø 26 cm, H 19,5 cm) allows you to easily cure multiple products at once. In addition, the box has enough space to place an articulator inside. If necessary, the UV light bulbs can be changed easily.

The new LC-3DPrint Box allows you to post-cure NextDent 3D materials. Post-curing is required in order to obtain the final material properties. After the printing process the printed material is not yet fully cured and should be treated. First, the material needs to be cleaned in an alcohol bath followed by the post-curing cycle. Post-curing is an UV light treatment to ensure that NextDent materials obtain full polymer conversion, through this the residual monomer is reduced to a minimum and the highest mechanical properties are obtained. This procedure is a necessary step to produce a biocompatible end-product. Always follow the instructions for use relevant to the corresponding material.

Specification	Value
Related voltage	110V/230V,50/60Hz,2,6A/1,3A
Powerconsumption	10 W
Fuse	T2.0A,AC250V
Dimension	W41xL44xH38cm
Weight	22 Kg



LC-3D Mixer

The LC-3D Mixer is a roller/tilting stirring device for mixing 3D printing materials before pouring in the Resin Tray of the printer. Print resins must be mixed well. Only handshaking is insufficient for highly filled and colored materials. The LC-3D Mixer should be used to thoroughly mix the material. When mixed insufficiently color deviation and print failures may occur. The NextDent LC-3D Mixer keeps your NextDent 3D materials ready for use at any time at an optimum consistency.

Specification	Value
Related voltage	AC 100-240V, 50/60 HZ
Powerconsumption	10 W
Fuse	250 V, T2A
Dimension	W410xD270xH100mm
Weight	4Kg

Overview NextDent Materials



Denture 3D+ 3D Printresinforthe manufacturing of removable denture bases



C&B MFHMicro Filled Hybrid 3D Print resin for the manufacturing of long-term temporaries



Try-In
3D Print resinfor the manufacturing of try-in devices



Tray3D Printresinfor
the manufacturing of
individual impression trays



SG 3D Printresin for the manufacturing of dental surgical guides



Model 2.0
3D Print Resinfor the manufacturing of prosthodontic and orthodontic models



Gingiva Mask
3D Printresinfor
the manufacturing of gingiva
masks on dental models



a Mask Model Ortho

resinfor 3D Printresinfor

uring of gingiva the manufacturing
ental models of orthodontic models



Ortho Clear
3D Printresinfor
the manufacturing of dental
splints and retainers



Ortho Rigid 3D Printresin for the manufacturing of dental splints



Ortho Indirect Bonding Tray 3D Printresinfor themanufacturing of Indirectbonding trays



Cast 3D Printresin for the manufacturing of castable parts



NextDent 5100

NexDent51003D printer powered by revolutionary Figure4™ technology



LC-3DPrint Box High-capacityUV postcuring unit



LC-3DMixer
Roller/Tilting mixer
for the mixing of
3D printing materials



by 🔱 3D SYSTEMS





U.A.E Office : 3204, Prism Tower, Business Bay, P.O. Box 28820, Dubai, U.A.E.

K.S.A Office: Al Saif Tower, 6th Floor, King Abdullah Street, Near Intersection with Olaya Street, Riyadh, K.S.A.

Tel:+971.4.443.3853; Fax:+971.4.443.3938; Email:info@3d-me.com; Website: www.3d-me.com